

MATERIALS LIST

Mat. ID	CE	Material Name	Brand	Ingredients
6	0123	IPS E-MAX CERAM	IVOCLAR	SiO ₂ , Al ₂ O ₃ , ZnO ₂ , Na ₂ O, K ₂ O, ZrO, CaO, P ₂ O ₅ , fluoride, pigments
8	0123	IPS E-MAX PRESS INGOTS	IVOCLAR	SiO ₂ , Li ₂ O, K ₂ O, MgO, ZnO, Al ₂ O ₃ , P ₂ O ₅
10	0120	D-sign 98	IVOCLAR	Au85.9%, Pt12.1%, Zn1.5%, In<1.0%, Ir<1.0%
12	0120	Lodestar	IVOCLAR	Au51.5%, Pd38.5%, In8.5%, Ga1.5%, Ru <1.0%, Re <1.0%
13	0120	Evolution Lite	IVOCLAR	Au40.3%, Pd39.3%, Ag9.2%, In9.3%, Ga1.8%, Ru<1.0%
14	0120	MiniGold	IVOCLAR	Au40.0%, Pd4.0%, Ag47.0%, Cu7.5%, Zn1.0%, In<1.0%, Ir<1.0%
17	0197	BEGO WIRONIT extrahart	BEGO	Co 63%, Cr 30%, Mo 5%, Si 1.1%, Mn, C
18	1275	Girobond [®] nbs	Girrbach	Co62.4%, Cr 25.5%, Mo 5.1%, W5.2%, Si1.1%, Nb, Fe, N<1%
19	0120	Self Curing-Pink	Vertex	Polymethylmethacrylate, plasticiser, pigments and catalysts, 0.5-1.5% Benzoylperoxide
20	0483	Orthocryl Powder-Clear	Dentaurum	Methylmethacrylate95%, Butandiole, Dimethacrylate<5%
22	0123	Ivostar	IVOCLAR	Methylmethacrylate
23	0483	Ropak uv-fluid	BREDENT	
25	0044	M.L. Primer	SHOFU	Azeton80%, other
26	0470	Preci-Vertex	CEKA	
27	0483	FINO AUTOBAOND	FINO	ATT Adhesive
28	0470	VALPLAST	VALPLAST	Nylon
29	0120	HIGH FUSING WHITE CERAMIC SOLDER	IVOCLAR	Au45.0%, Pd12.4%, Ag41.5%, In1.0%, Li<1.0%, Ru<1.0%
30	0120	SPECIAL HIGH FUSING WHITE CERAMIC SOLDER	IVOCLAR	Au47.0%, Pd10.3%, Ag41.0%, B<1.0%, In1.4%, Ru<1.0%, Other<1.0%
31	0120	UNIVERSAL SOLER 1015W	IVOCLAR	Au18.5%, Pd5.99%, Ag72.5%, In3.0%, Ir<1.0%
32	0120	UNIVERSAL SOLDER PKF	IVOCLAR	Au48.8%, Ag40.5%, Zn7.3%, Pt2.8%, Other<1.0%
38	0123	E-MAX GLAZE	IVOCLAR	Ethandiol
39	0123	E-MAX CERAM MARGIN	IVOCLAR	SiO ₂ , Al ₂ O ₃ , ZnO ₂ , Na ₂ O, K ₂ O, ZrO, CaO, P ₂ O ₅ , fluoride, pigments

Mat. ID	CE	Material Name	Brand	Ingredients
53	0483	MK1	DENTAL-ATTACHMEN	
55	0044	TK-SNAP	SI-TEC	
56	0123	E-MAX CERAM GINGIVA	IVOCLAR	SiO ₂ ,Al ₂ O ₃ ,ZnO ₂ ,Na ₂ O,K ₂ O,ZrO, CaO,P ₂ O ₅ ,fluoride,pigments
58	0124	TITANIUM CAD CAM ABUTMENT	GEO	Titanium80-100%,Molybdenum1-15%, Aluminum<10%, Zirconium<10%,Iron<5%,Chromium<5%,Copper<5%,Tin<5%,Cobalt0-Trace,Manganese0-Trace, Nickel0-Trace,Tantalum0-Trace, Tungsten0-Trace,Vanadium0-10%
60	1251	BALL ATTACHMENT	CEKA	
63	0124	TITANIUM CERAMIC	VITA	SiO ₂ 63-68%,Na ₂ O 5-7%,Cao 1.5-2.5%,K ₂ O 8-10%,Al ₂ O ₃ 8-10%,BaO 0.3-1%,ZrO ₂ <0.5%,SnO ₂ <0.5%,B ₂ O ₃ 7-9%
64	0120	PMMA DISK	YAMAHACHI	
66	0483	DENTAURUM CLASP	DENTAURUM	NICKEL FREE
71	0483	BALL RETAINER CLASPS	DENTAURUM	
72	0197	Upzir solid zirconia	UPCERA	ZrO ₂ +HfO ₂ +Y ₂ O ₃ ≥ 99%,Y ₂ O ₃ =4.5~6.0%,HfO ₂ ≤ 5%,Al ₂ O ₃ ≤ 0.5%,other ≤ 0.5%
74		Resilient Acrylic	Bosworth Trusoft	
77	0297	mega TRIM Provi Implant	megadental	
78	0483	Remanium Draht/Wire	Dentaurum	
79	0483	Expansion Screw	Dentaurum	
80	0344	CZR	NORITAKE	
81	0120	Harmony Hard	IVOCLAR	Au74.0%,Pd3.8%,Ag12.0%,Cu9.0%,Zn<1.0%, In<1.0%, Ir<1.0%
82	0120	SPARTAN PLUS	IVOCLAR	Au2.0, Pd78.8, Cu10.0, Ga9.0, Ir<1.0, Ge<1.0, Li<1.0
83	93/42/EWG	Silensor-sl Teilekarten	ERKODENT	polyethylene
85	93/42/EWG	ERKODUR	ERKODENT	Ethylene terephthalate,Ethylene glycol,Vinyl acetate, Ethylene polymerization
86	93/42/EWG	ERKOFLEX-BLEACH	ERKODENT	Ethylene terephthalate,Ethylene glycol,Vinyl acetate, Ethylene polymerization
87	0482	DDcubeX2	Dental Direkt materials	ZrO ₂ +HfO ₂ +Y ₂ O ₃ ≥ 99%; Y ₂ O ₃ < 10%; Al ₂ O ₃ ~ 0.05%; other ≤ 0.05%
88	289250	Vario Ball Snap	BREDENT	
90	0123	IPS Classic	IVOCLAR	SiO ₂ 40~65wt.%, Al ₂ O ₃ ,B ₂ O ₃ , BaO, CaO, CeO ₂ , K ₂ O,MgO, Na ₂ O, P ₂ O ₅ , SnO ₂ , TiO ₂ , ZrO ₂ and pigments
91	0123	SR Nexco Opaquer	IVOCLAR	IVOCLAR VIVADENT

Mat. ID	CE	Material Name	Brand	Ingredients
92	0123	SR Nexco Paste	IVOCLAR	IVOCLAR VIVADENT
95	0123	E-MAX GLAZE	IVOCLAR	Ethandiol
96	0123	IPS E-MAX CERAM GLAZE SPRAY	IVOCLAR	Ethandiol
97	0124	Reinforcement mesh nets	RENFERT	
98	0120	THERMOSENS RIGID	VERTEX	Thermoplastic Dental Resin
102	524151MR5	GIALLOY TITANIUM ALLOY BLANKS	SRL	Ti 90.04%, Al5.94%, V4.0%,Fe0.02%,O 0%