

Compatible to Implant Systems	Platform	Ti-base	Premill	Analog
Osstem Implant TS	TS Mini 3.0/3.5	✓	✓	✓
	TS Standart 4.0/4.5/5.0	✓	✓	✓
Osstem Implant SS	SS Regular 4.8	✓	✓	✓
Dentium Super-Line Implantium	Universal 3.7/4.0/4.5	✓	✓	✓
Megagen AnyRidge	Universal 3.5/4.0/4.5	✓	✓	✓
Megagen AnyOne	Universal 3.0/4.0/4.5/5.0	✓	✓	✓
Dio UF / UFII	Narrow 3.0/3.3		✓	
	Standart 3.8/4.0/4.5/5.0	✓	✓	✓
NeoBiotech IS	Standart 3.5/4.0/4.5/5.0	✓	✓	✓
Straumann Bone Level	SC 2.9			
	NC 3.3	✓	✓	✓
	RC 4.1/4.8	✓	✓	✓
Straumann Tissue Level	NN 3.5			
	RN 4.8	✓	✓	✓
	WN 6.5	✓	✓	✓
Nobel Biocare Active	3.0			
	NP 3.5	✓	✓	✓
	RP 4.3/5.0	✓	✓	✓
	WP 5.5	✓	✓	✓
Nobel Biocare Branemark	NP 3.5	✓	✓	✓
	RP 4.1	✓	✓	✓
	WP 5.1	✓	✓	✓
Noble Biocare Replace	NP 3.5	✓	✓	✓
	RP 4.3	✓	✓	✓
	WP 5.0	✓	✓	✓
	6.0	✓	✓	✓
Dentsply Xive/Frialit 2	3.0	✓	✓	✓
	3.4	✓	✓	✓
	3.8	✓	✓	✓
	4.5	✓	✓	✓
	5.5	✓	✓	✓

Compatible to Implant Systems	Platform	Ti-base	Premill	Analog
Dentsply Ankylos	Universal 3.4/4.5/5.5	✓	✓	✓
Biomet 3i Certain	3.4	✓	✓	✓
	4.1	✓	✓	✓
	5.0	✓	✓	✓
	6.0	✓	✓	✓
Biomet 3i Osseotite	3.4	✓	✓	
	4.1	✓	✓	✓
	5.0	✓	✓	
	6.0	✓	✓	
Implant Direct 3.0	3.0	✓	✓	✓
Zimmer - Tapered Screw-Vent BioHorizons - Tapered Internal Implant Direct - Legacy BlueSkyBio - Internal	3.5	✓	✓	✓
	4.5	✓	✓	✓
	5.7	✓	✓	✓
Astra Tech Osseospeed	X-Small 3.0	✓	✓	
	Aqua 3.5/4.0	✓	✓	✓
	Lilac 4.5/5.0	✓	✓	✓
Astra EV System	3.0	✓	✓	✓
	3.6	✓	✓	✓
	4.2	✓	✓	✓
	4.8	✓	✓	✓
	5.4	✓	✓	✓
Camlog	3.3	✓	✓	✓
	3.8	✓	✓	✓
	4.3	✓	✓	✓
	5.0	✓	✓	
	6.0	✓	✓	
MIS SEVEN	Narrow 3.3	✓	✓	✓
	Standart 3.75/4.2			
	Wide 5.0/6.0	✓	✓	✓
MIS C1	3.3	✓	✓	✓
	3.75/4.2	✓	✓	✓
	5.0		✓	✓

Compatible to Implant Systems	Platform	Ti-base	Premill	Analog
Thommen SPI	3.5			
	4.0	✓	✓	✓
	4.5	✓	✓	✓
	5.0			
	6.0			
Biohorizons Laser-Lok 3.0	3.0	✓	✓	✓
Anthogyr AXIOM	Universal 3.4/4.0/4.6/5.2	✓	✓	✓
Neodent GM	Standart 3.5/3.75/4.0/4.3/5.0		✓	
Bicon Implant	2.0			
	2.5			
	3.0	✓	✓	
Medentis medical ICX-System	Universal 3.75/4.1/4.8	✓	✓	